

12D MAX PLUS



Innovative technology rejuvenate youth
Bidirectional output



With our unique
patented technology
the results are



Brighten
skin tone



Delicate skin



Regenerated
collagen



Smooth
fine lines



Enhance
elasticity



Fat pad
reduction



Improved sexual
happiness index




Improve
firmness



ULTRA MICRO PULSE ULTRASONIC SYSTEM

Hall high-frequency brushless motor works in tandem (left and right) to ensure efficient operation and smooth operation. The high rate CPU processor makes the output waveform smoother and more unique. The voltage can belowered to drive the high-frequency energy output, so that the focusing energy is more stable, more concentrated and more efficient, so that you can enjoy safely. Dual AI interfaces and dual gamepads work together to achieve the effect that 1+1 is greater than 2.



INNOVATIVE TECHNOLOGY


Two-way output technique Active efficient operation and precise circuit operation	MP Pattern (Multi-point line) In MP mode, the microsecond energy transmission mode has smaller and denser gathering points, higher security points are distributed more evenly, and the effectiveness is stronger	Large focal areas gather energy Effectively reduce the ultrasonic energy density in the treatment area, make the heat effect more uniform and gentle, and greatly reduce the side effects	Multi-mode combined nursing High frequency ultrasound treatment on the skin layer to the point, surface, line and area layer
--	--	--	---

Superpulse ultrasound | In a very short period of time, the release of extremely high pulse energy, the action rate is several times that of the ordinary treatment mode

TECHNOLOGY

HIGH FREQUENCY FORWARD AND REVERSE LINKAGE WORKDUAL HANDLE COLLABORATIVE WORK


In MP mode, the microsecond energy transmission mode has smaller and denser gathering points, higher security points are distributed more evenly, and the effectiveness is stronger



Nomal Model nursing

Two-way output technique

17 points / 1 engine in total



MP Model nursing

Two-way output technique


417 points / 1 engine in total



Aggregate energy effect

Two-way output technique

Large focal areas gather energy



Multiple frequency probe

Multi-mode combined nursing




Nursing time

Superpulse ultrasound

Mechanical effect + thermal effect + cavitation effect

Using ultra-micro pulse ultrasonic system, multi-frequency ultrasonic treatment probe on the point, surface, line and skin layer by layer! Maximum guarantee of curative effect. The concentrated energy in the large focal area effectively reduces the ultrasonic energy density in the treatment area, making the heat effect more uniform and gentle, and greatly improving the experience. Energy releases very high peak energy in a very short period of time, and the action time is shorter, and the treatment time is shortened.



TREAT VAGINAL RELAXATION AND RESTORE TIGHTNESS, ELASTICITY AND SMOOTHNESS



TIGHT

Tighten vagina
Treat loose vagina

ELASTIC

Revitalizes skin
and enhances elasticity

SMOOTH

Sensitive upgrade to enjoy
moisturizing

LEFT Superpulse handle

The focal probe is applied to the skin

- 1.5** Face - 1.5 - 5.0MHz
For forehead, nose, chin, lips, cheeks, skin tone, skin pores
- 3.0** Face - 3.0 - 5.0MHz
To treat the forehead and stimulate collagen growth
- 4.5** Face - 4.5 - 5.0MHz
Focus the focal point of 0.5MM and reduce the muscle layer to achieve slimming effect
- 6.0** Cheek - 6.0 - 5.0MHz
For the cheekbone, skin tone, skin pores and pore shrinkage



RIGHT

RIGHT Superpulse handle

The body probe acts on the skin

Use frequency treatment method on the skin layer to be treated. Focus the focal point of the treatment area to achieve the desired effect.

- 1.5** Eye - 1.5 - 5.0MHz
Treating orbital cellulose and problems such as double eyelids, skin tone, and eye bags
- 3.0** Body - 3.0 - 5.0MHz
Treat the whole body, skin tone, skin pores, and skin texture
- 12.0** Body - 12.0 - 5.0MHz
To improve the body and stimulate collagen growth
- 18.0** Body - 18.0 - 5.0MHz
Focus the focal point of 1.5MM and reduce the muscle layer to achieve a slimming effect

PRIVATE HEAD 3.0MM

The energy reaches the submucosal tissue directly, accelerating cell activity and causing collagen regeneration and reorganization, increasing mucosal elasticity and tightening vaginal muscles.

PRIVATE HEAD 4.5MM

The energy reaches the fascia layer directly, causing the fascia layer to thermally solidify and produce the effect of tightening the vagina, improving the muscle structure supporting the mucosa from deep to shallow, and better helping the muscle layer to restore elasticity and become tighter.



PRIVATE TEST TUBE

Adapt the airbag pressure measurement principle to detect vaginal relaxation



PRIVATE CUSTOM FULL LAYER ANTI-AGING

According to the age group, skin status, customer needs, point, line, surface combination, to achieve a full-layer, all-round anti-aging private customized solution.



- Epidermal and superficial dermal activation.
- Deep dermal activation.
- Tighten the fascia instantly.
- Enhance the supporting force of fascia layer and reduce the slack tissue.
- Enhance cell activity and increase cell water storage capacity.
- Increase facial volume, stimulate collagen regeneration, elastic fiber network regeneration, full, elastic.




SINGLE MODE NURSING

MP mode (multi-point line formation)
[417 points/1 engine in total]

Dot mode can also be selected
[17 points/1 round in total]

The mode can be changed freely by touching the interface of the same probe



MULTI-MODE COMBINED NURSING

The micro-pulse realizes the superposition of tiny and precise thermal condensation points, and uniformly and intensively acts on all areas of the target level. The dense and scattered energy acts more widely, highlighting the overall advantages.

HIGH-INTENSITY FOCUSED ULTRASOUND

High-energy focused ultrasound reaches the vaginal SMAS layer directly, lifting the SMAS fascia suspension and completely solving the problem of vaginal relaxation



13.3-INCH HD SMART SCREEN



High frequency forward and reverse linkage work
Dual Handle Collaborative Work



VISUAL OPERATION INTERFACE

No blind pressing, visual operation interface, more intuitive to use





ARGUMENT

PRODUCT CONFIGURATION

PRODUCT NAME	12D-MAX plus
MODEL	12D-MAX
POWER SUPPLY	AC 110V AC 220V
POWER	650W
FREQUENCY	2 MHz 4 MHz 5.5 MHz
INSTRUMENT WEIGHT	10 kg
GROSS PRODUCT QUANTITY	20 kg
PACKAGE SIZE	60 X 43 X 57 cm



PRECAUTIONS AFTER OPERATION

- 1 It is normal to have temporary slight redness and swelling after operation; Pay attention not to go to high temperature environment such as hot sauna yoga and strenuous exercise within 1 week after operation, and do not expose to the sun;
- 2 Try not to make up on the day after operation (if makeup needs to be removed, pay attention to the strength, and gently remove makeup oil to clean the face);
- 3 Do not clean your face with overheated water within a week (water that does not exceed body temperature is enough); Pay attention to water supplement. It is recommended to apply water supplement facial mask every morning and evening for 3 days or one week; Skin care and moisturizing products should be wiped 2-3 times as much as usual;
- 4 Don't eat anything with heavy color such as soy sauce within a week, and forbid alcohol, spicy, seafood and other foods;
- 5 It is not recommended to operate similar instruments at the same place on the same day after operation;
- 6 After operation, it is recommended to use laser and other instruments after skin repair. (The same part shall be adjusted according to the skin condition of the day);
- 7 There may also be a small white energy rash or skin erythema at some operating sites, which are normal, mild and temporary reactions (can be treated with ice or growth factor repair products according to the specific situation)

ZENLASE

Aesthetic Laser Equipment · Factory Direct

WHATSAPP

+86 151 0211 7919

EMAIL

info@zenlase.com

WEB

www.zenlase.com

HOURS

Mon-Sat 09:00-18:00 (GMT+8)

Quotations, orders, spare parts and after-sales for overseas buyers are handled by the ZENLASE team.