



MEDICAL AESTHETIC SOLUTIONS

QM - 4S





THE COOLEST WAY TO LOSE WEIGHT

■ What will my results be?

Each Cryolipolysis treatment results in a significant noticeable reduction of fat in the treated area, typically between 20% to 40%.

■ When will I see results?

You may start to see changes as quickly as three weeks after your treatment, and will experience the most dramatic results after two months. But your body will still flush out fat cells and continues doing so for up to four to six months after treatment.

SEE A SLIMMER YOU



Technical Specifications

Input voltage	AC110V/220V 50-60Hz
Power consumption	< 500w
Screen	10-inch wide color touch screen
Cooling device output temperatur	-11-45
Cooling device output pressure	0-100KPa
Cooling liquid	pure water
Default time	0-120min
Fuse size	F2AL250V
Ambient temperature	5 -+40
Relative humidity	80%
Atmospheric pressure	80KPa-106KPa

GIVE YOU MORE DIFFERENT DEMONSTRATION

✓ Four Different Size Cryolipolysis Handle

Four different sizes of handles which can apply to different parts of body according to the fat thickness of each client.

✓ Two Applicators Working Simultaneously

About our 100 size, 150 size, 200 size cryolipolysis treatment handles. Any 2 handles (except the 300 size cryo handle) can work at the same time. Save more of your treatment time.

✓ Touch Screen on Each Applicator

Touch screen on the applicator makes it easy to adjust the setting.

✓ Double Air Pumps

When two applicators working simultaneously, a single pump is not enough, our machine equipped with two separate powerful air pumps to guarantee enough vacuum pressure.

✓ High Quality Material & Nice Workmanship

The machine and applicators are made of high quality ABS. Our assemblers are well-trained and responsible for the quality of whatever they make.

✓ Beautiful Exterior Design

With nice silver streamline housing and great details, it's sure the machine will help improve the image of your salons or clinics



Longer cooling plates

Maximizes treatment area with 38% greater cooling surface

Flat applicator cup



Push-button switch control

Patient-friendly design



Unique applicator design



Touch screen

Touch screen on the applicator makes it easy to adjust the setting.

Thickness	Dose	Time
<20mm	0 °C	60min
20-40mm	-5 °C	45min
30-50mm	-8 °C	40min
50-60mm	-10 °C	40min
>60mm	-11 °C	40min

Freeze -11°C

Strong semiconductor cooling system



Red LED light

Speed up the blood circulation and oxygen to increase activity



Green LED light

Reduce the skin swelling in the freezing treatment

	Model	Size
Criolipolysis treatment handle	Size 100	90*43*65mm
	Size 150	130*50*65mm
	Size 200	190*70*92mm
	Size 300	350*120*120mm
Antifreeze membrane	ETG I	280*280mm
	ETG II	270*300mm
	ETG III	340*420mm



■ **SIZE 100**

- Abdomen and flanks, also suggested for the arms and inner thighs



■ **SIZE 150**

- Fit on curved or narrow parts of the body



■ **SIZE 200**

- Best matches with the contour of the abdomen




■ **SIZE 300**

- Suggested for large volume reduction or de-bulking



WHAT CAN CRYOLIPOLYSIS REDEFINE ON YOU?

Eliminates 20% to 40% of Your Fat with ONE Single Treatment



ARE YOU SATISFIED WITH YOUR ANTIFREEZE MEMBRANE?



✓ **SIZE: 34*42CM WEIGHT: 110G**

✓ **SIZE: 27*30CM WEIGHT: 70G**

✓ **SIZE: 28*28CM WEIGHT: 60G**



ZENLASE

Aesthetic Laser Equipment · Factory Direct

WHATSAPP

+86 151 0211 7919

EMAIL

info@zenlase.com

WEB

www.zenlase.com

HOURS

Mon-Sat 09:00-18:00 (GMT+8)

Quotations, orders, spare parts and after-sales for overseas buyers are handled by the ZENLASE team.