



Innovative portable IPL SHR System!

Model: QM-100+



Applications:

- Hair removal
- Skin rejuvenation
- Pigments removal
- Acnes Removal
- Vascular Removal
- Wrinkle removal




Working Principle:

It is a revolutionary approach, rather than bombarding your skin with high amounts of damaging energy. SHR fires multiple shots but at low Joules, in doing so it gently heats the hair follicle to the heat needed and the most you will feel is warmth and a tingling sensation, some clients comparing it to a warm massage. SHR also utilizes In-Motion technology, where the hand piece is always in motion over the skin.



Specialty:

© Smart and compact metal design, it can separate into 2 parts, which solve the Customs questions in some countries., and reduce the complexity of after-sale service.



© Professional Ergonomic design Handpiece.

© Innovative filters for different treatments



Unique and Elegant Design:



Technical Specifications

Item	IPL	SHR
Model	QM-100+	QM-100+
IPL Peak Power	1500W	2000W
Wavelength (Spectrum)	430/530/640—1200nm	650-950nm
Energy Density (Fluence)	10-60J/cm ²	1-26J/cm ²
Spot Size	15*50mm	15*50mm
Pulse Repetition Rate	0.5-2Hz	1-10Hz
Pulse Duration	1-50ms	1-15ms
Cooling	+Continuous Crystal contact cooling (-5℃~1℃) +Air cooling +Closed water circulation cooling	+Continuous Crystal contact cooling (-5℃~1℃) +Air cooling +Closed water circulation cooling
Pulses	1-6	Single Pulse
Stand-by Working	20 hours	20 hours
Display	10.4" True Color LCD Touch Screen	
Electrical Requirements	110/230VAC, 15/20A max., 50/60Hz	

Before & After



ZENLASE

Aesthetic Laser Equipment · Factory Direct

WHATSAPP

+86 151 0211 7919

EMAIL

info@zenlase.com

WEB

www.zenlase.com

HOURS

Mon-Sat 09:00-18:00 (GMT+8)

Quotations, orders, spare parts and after-sales for overseas buyers are handled by the ZENLASE team.