

HIFU-Q3

Let Miracles Happen !



Skin rejuvenation

Skin tightening and lifting

Wrinkle removal

Improve skin laxity

Recover skin elasticity


Applications

- Skin rejuvenation
- Skin tightening and lifting
- Wrinkle removal
- Improve skin laxity
- Recover skin elasticity




Working Principles


HIFU is a non surgical and a non laser procedure that directs high intensity focused ultrasound also known as hifu knife into the deep skin layer in order to gently heat the tissue to lift and tighten aged skin without any incisions or disrupting skin surface .Hifu technology is able to transfer the heat energy into 3-4.5mm under the surface of the skin .It will heat the smas up to 60-70 temperature the raises thermal coagulation point which accelerates skin rejuvenation process and maximizes tightening effect of the skin layer .



Advantages of HIFU-Q3

1. Adopts HIFU technology, the temperature of focused ultrasound can reach 65 °C to 100 °C
2. The treatment by HIFU is non-invasive and painless, HIFU-Q3 can works at areas that 3mm and 4.5mm underneath of skin surface, focus at SMAS, stimulates collagens increase, recover SMAS's function of stretch the muscles and skin, get a perfect result of freckle removal.
3. HIFU-Q3 has extra transducer 13mm give you an option for body treatment: body tightening ,lifting and slimming.
4. High accuracy and controllability, with the unique ultrasound skin scan system, HIFU-Q3 can lock the target area for treating precisely, easy for targeted therapy, parameters are adjustable for safe, accurate treatment
4. Non-invasive treatment, without any downtime .Because strong energy is delivered deeply into skin without affecting the epidermis the treatment is safe .
5. Minimization of treatment affects because facial cleansing and makeup is possible even immediately after treatment .





Comprehensive treatment





ZENLASE

Aesthetic Laser Equipment · Factory Direct

WHATSAPP

+86 151 0211 7919

EMAIL

info@zenlase.com

WEB

www.zenlase.com

HOURS

Mon-Sat 09:00-18:00 (GMT+8)

Quotations, orders, spare parts and after-sales for overseas buyers are handled by the ZENLASE team.