

— Smart Tecar Wave —



Shock Wave



EMS



Tecar



Electromagnetic shockwave

SHOCKWAVE is aimed at the affected areas that are the source of chronic pain. The influence of the shockwaves causes to the dissolution of calcium deposits and leads to better vascularization. The after-effect is relief from the pain.



Handpiece:


Tips:



1#



2#



3#



4#



5#



6#



7#

QMmed
MEDICAL AESTHETIC SOLUTIONS

Beijing QM Medical Co., Ltd.
www.zenlase.com

Tecar :Capactive And Resistive Energy Transfer



TECAR therapy as a system of Capacitive and Resistive electric transfer ,which is one of the methods used in diathermy, was developed as a form of deep thermotherapy ,delivers radiofrequency (RF) energy,which passes between active electrode and inactive electrode, and generates heat in the human body .

Handpiece:

Optional
80MM
Handpiece:


Tips:



CET Handle



RET Handle



Large/Middle/Small Sizes

Electric Muscle Stimulation


EMS “Electrical Muscle Stimulation” is an internationally accepted and proven way of treating muscular injuries. It works by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively. When the muscle receives this signal, it contracts as if the brain has sent the signal itself. As the signal strength increases, the muscle flexes as in physical exercise. Then when the pulse ceases, the muscle relaxes and the cycle starts over again, (Stimulation, Contraction and Relaxation.)



EMS Cup:



Main Interface



EMS SHOCK Interface



Tecar Interface

Resistive Energy Transfer


RET technique works in tissues with higher resistance such as bones, tendons, and joints. The deeper penetration of RET also helps dissolve and soften deep stubborn fat and circulates the fat out of the body to achieve weight loss.



Capasitive Energy Transfer

CET technique works in tissue containing a high content of electrolytes such as muscle and soft tissue.

Smart Tecar Wave



Feature



Physical Therapy

Orthopaedic disorders of the peripheral nervous system, Pelvic floor rehabilitation
Acute and chronic pain...

Rehabilitation

To speed up muscle recovery and help athletes return to play as soon as possible...

Pain Management

Treats injuries of knee, shoulder, hip, ankle, spine, hands and muscles...

Decrease abdominal and body fat


ED treatment (Shock wave)

Muscle growth strength training post-exercise recovery (EMS)




Shockwave for physiotherapy

- Shockwave is an acoustic wave which carries high energy to painful areas and soft tissues with subacute, subchronic and chronic conditions. This energy promotes healing, regenerating and reparative processes. It's a unique, non-invasive solution for pain associated with the musculoskeletal system. Just three to four treatments needed at weekly intervals. A therapy session only takes approximately 10 minutes.



ED Treatment

Most men suffering from erectile dysfunction have vascular problems affecting the vessels that supply blood to the cavernous bodies of the penis, resulting in a decreased ability to develop and maintain an erection. The end result is an inadequate blood-flow followed by premature draining of an already reduced amount of blood from the cavernous bodies, meaning an erection is not achieved. Shockwave Therapy for ED of this type can be a highly effective treatment.



Tecar Effect

The heat speeds up metabolism. This causes blood to flow faster and to become more oxygenated. The result is that more oxygen, and other healing properties of your body's natural systems, are rushed to the site. The waste is also removed more quickly. The overall result is that your pain is reduced significantly, and the injury is healed more quickly.




Muscle growth strength training post-exercise recovery



EMS (electrical muscle stimulation) is a machine that delivers a stimulating pulse to your muscles.

EMS builds and grows muscles by causing this contraction. When the muscle is tightened it is performing work, and this causes lactic acid to eventually enter the muscle, this can lead to faster muscle growth.



QMmed

MEDICAL AESTHETIC SOLUTIONS

| Beijing OM Medical Co., Ltd.

www.zenlase.com

ZENLASE

Aesthetic Laser Equipment · Factory Direct

WHATSAPP

+86 151 0211 7919

EMAIL

info@zenlase.com

WEB

www.zenlase.com

HOURS

Mon-Sat 09:00-18:00 (GMT+8)

Quotations, orders, spare parts and after-sales for overseas buyers are handled by the ZENLASE team.